



Watch the instructional video:
<https://alere.wistia.com/medias/guym62w5j>

Revision Date: 2022-03-31

INSTRUCTIONS FOR USE

Catalog Number: 29012-W01

EN

Before testing you must read all the steps. Conformance with the test procedure is necessary to ensure an accurate result.

Precautions


Do not use

- If you have a bleeding disorder
- If you are on HIV treatment (ARVs)
- If you are needle phobic
- If the kit bag or components is broken
- If the kit or components have been used
- If your area is under poor lighting

Do not eat or drink while you perform the test


STEP 1: PREPARATION

1 Prepare a Timer and Tissue.



⚠ Not included in the kit, but needed.

2 Open and place all materials on a flat and clean surface with bright light.



3 Wash hands in warm water and dry. If no warm water is available, rub your hands together.




4 Choose ring finger or middle finger.




Avoid dominant hand.

STEP 2: COLLECT BLOOD


5 Massage and rub your hand & finger to increase circulation.




6 Clean your finger with B Alcohol Pad. Let it dry for 10 seconds.




7 Remove the Lancet cover.




8 Press the F Lancet against the finger until it clicks.



9 Massage from the base to the tip, let 2 drops of blood fall into the Basin. If you are having difficulty, wipe finger clean and squeeze again.




10 Wipe finger with B Alcohol Pad and apply the C Plaster. If needed, press on the plaster to stop bleeding. Start next step immediately to transfer blood.




STEP 3: TEST


11 Squeeze the top of the E Specimen Dropper all the way down and hold while dipping into the blood sample.



Dip the dropper into the blood in Basin and release slowly to draw blood into the dropper. ⚠ Avoid bubbles when drawing blood.




Place dropper over the Well labelled S. Squeeze the top of the dropper to apply 1 drop of blood into the Well labelled S.



12 Hold D Buffer bottle vertically and apply 1 drop of Buffer into the Well labelled S.



13 Start the Timer. Read the result in 15-20 minutes, do not read past 20 minutes.




See flip side for STEP 4:
Read Result

Number: 1156182105



STEP 4: READ RESULT

14 Read test result in 15-20 minutes.



INTENDED USE

The CheckNOW™ HIV SELF TEST is a single-use, *in vitro* (outside the body) visually read rapid immunoassay that uses a blood sample from a finger puncture for the qualitative detection of antibodies to HIV-1 and HIV-2 in blood. The CheckNOW™ HIV SELF TEST is intended to be used manually by untrained lay users (self testing) who are 14 years and older to aid in the diagnosis of HIV-1 or HIV-2 infection. This test is not intended to be used as an HIV screening test for blood donation.

SUMMARY

HIV is recognized as the virus that causes AIDS (Acquired Immunodeficiency Syndrome). The virus is transmitted by sexual contact, exposure to infected blood, certain body fluids or tissues, and from mother to child during pregnancy¹. The CheckNOW™ HIV SELF TEST detects the presence of antibodies to HIV-1 and/or HIV-2 in blood. The product includes a Test device and a Buffer. To use the test, two drops of blood sample are collected from the fingerstick in the Basin of the plastic cover. One drop of blood is transferred by a Specimen dropper to the Well. After that, one drop of Buffer is applied. When the test is completed, two lines can appear on the device. The red line in the Control Line (C) area will only become visible if the added blood sample and/or buffer have moved over the T/C Line areas of the reading window. The T line area is pre-coated with HIV-1 antigen glycoprotein 41 and HIV-2 antigen glycoprotein 36. The red line in the Test Line (T) area will only become visible if the applied sample contains antibodies to HIV-1 or HIV-2.

STORAGE

- Store the test kit at 2-30°C (36-86° F) until the expiry date. Do not freeze.
- Do not use the test kit after expiry date printed on the test pouch.
- Do not open the sealed foil test pouch until you are ready to use the test.
- Test device should be used within one hour after the pouch has been opened.
- The buffer is for single-use, and should be used within one hour after cap open.

WARNINGS AND PRECAUTIONS

- This test may give an unexpected positive result. Whether the result is positive or negative, you should consult with your doctor before making medical decisions.
- For *in vitro* diagnostic use only.

30°C	Store between 2-30°C		Consult instructions for use		Do not reuse		Manufacturer
	Catalogue number		Use-by date		Batch code		Sterilized using irradiation
	In vitro diagnostic medical device		Contains sufficient for <n> tests		Caution		Sterilized using ethylene oxide

NON-REACTIVE
(=Negative)



A line appears only in the C area. There is no line in the T area. The test did not detect the presence of HIV, however very recent exposure cannot be excluded.

It is recommended to conduct a retest after 6 weeks from latest risk of exposure to HIV.


REACTIVE
(=Positive)



One line in the C area, together with one line in the T area, no matter how faint, indicates the potential to be HIV positive.

Consult a health care provider. A reactive result must be confirmed by a lab test. Protect yourself and others! Avoid any activity that could transmit HIV to others.

TEST DID NOT WORK
(=Invalid)



If no line appears in the C area, even if a line appears in the T area, the test did not work.

Test again using a new kit or consult a health care provider.

- Perform test only by using a fresh blood sample. Not to be used with serum or plasma sample.
- The test device is recommended to be used at room temperature (15-30°C).
- Lancet and Specimen dropper are for single-use only. Do not reuse.
- If you do not understand the IFU please reach out to provided Technical Service contact.
- If your test is working, you will see a line in the Control Line area on your test device. If there is no line in the Control Line area, your test did not work, and the test result is invalid. However, the presence of a line in the Control Line area does not confirm sufficient specimen addition. You must transfer exactly 1 drop of blood to the Well before adding Buffer.
- Keep the SELF TEST and the components out of the reach of children.
- The buffer contains 0.09% sodium azide as a preservative, which may be toxic if ingested. If you get contaminated to your eyes rinse with running water for at least 15 minutes. If irritation persists, get medical attention.
- If your finger is still bleeding, use tissue or wipes. Protect yourself and others.
- Finger prick blood should be transferred from Basin to Well immediately to avoid blood clotting (within 2 minutes).

LIMITATIONS OF THE TEST

- The CheckNOW™ HIV SELF TEST is designed to be used with a finger puncture blood sample. Other body fluids must not be used.
- Not suitable for testing infant younger than 18 months².
- A NON-REACTIVE (Negative) result does not absolutely rule out the possibility of HIV infection.
- A REACTIVE (Positive) result must be confirmed by a health care provider using appropriate confirmatory testing.
- The intensity of the Test Line for a REACTIVE (Positive) result does not reflect how much HIV antibody is present in the blood sample.
- Although it's rare, false results may occur. If you have concerns that your result may be false, please contact your health care provider.
- Biotin concentrations up to 1500 ng/mL may lead to decreased Control Line intensity but have no impact on the internal control performance.

An incorrect or "false" NON-REACTIVE (Negative) result can occur for any of the following reasons: Incorrectly reading test result; Not following the Instructions for Use carefully; If you are on HIV treatment (ARV)^{3,4,5}; If you were very recently infected; The presence of bubbles during sample application, in particular in low positive samples.

An incorrect or "false" REACTIVE (Positive) result can occur for any of the following reasons: Incorrectly reading test result; Not following the Instructions for Use carefully; Having received an HIV vaccine; In cases of infection with cytomegalovirus.

TEST PERFORMANCE

The test has been shown in clinical evaluations performed by professional health care persons to correctly identify 99.9% (2097 out of 2100) with a confidence interval of 99.6% to 100% of HIV negative samples (known as the test's specificity). Further in field clinical evaluations conducted in South Africa, Congo, Vietnam and Spain, the test correctly identified 99.6% (1824 out of 1831) with a confidence interval of 99.2% to 99.9% of HIV negative samples when performed by first time self test users.

The test has also been shown in clinical evaluations performed by professional health care persons to correctly identify 100% (600 out of 600) with a confidence interval of 99.5% to 100% of HIV positive samples (known as the test's sensitivity). Further in field clinical evaluations conducted in South Africa, Congo, Vietnam and Spain, the test correctly identified 95.1% (270* out of 284) with a confidence interval of 91.9% to 97.3% of HIV positive samples when performed by first time self test users.

*Note: A total of 6 first time CheckNOW™ HIV SELF TEST users needed to be excluded from this analysis as they were observed to deny an unexpected result. Refer also to warnings and precautions about the use of this test.


To ensure that other medical conditions (potential cross-reactants) do not affect the performance of the CheckNOW™ HIV SELF TEST, samples of HIV negative blood were tested from people who had other conditions. These included 250 specimens from pregnant women and 342 other specimens as follows:

STEP 5: DISPOSAL

15 Place all used components back into the kit bag.



16 Seal the kit bag tightly.



17 Throw away the kit bag in waste bin or household rubbish.



Blood can transmit infectious diseases
Clean up spills

Dispose in accordance to local regulations

HAMA; Multiparous woman; Elevated IgG; Elevated IgM; Systemic lupus erythematosus; Hemolytic; Lipemic; Icteric; Rheumatoid Factor; ANA; Anti-E. coli positive specimens; Sickle-cell disease specimens; Blood from recipients of multiple blood transfusions; HBsAg; EBV; CMV; Malaria; Measles; Tuberculosis; Varicella zoster virus; Influenza A and B; Tick borne encephalitis; Influenza vaccine recipient; Human African trypanosomiasis; Yellow fever virus; Post-immunization measles; Vaccine-induced HIV seropositivity; Yellow fever vaccine recipient, Leishmaniasis positive; Syphilis; Toxoplasmosis; Helicobacter pylori; HSV; anti-HCV, anti-HBs, anti-HBc; anti-HTLV-1/2; anti-HEV, anti-HAV. These non-HIV medical conditions did not affect the performance of CheckNOW™ HIV SELF TEST with exception of the observed cross-reactivity seen with 2 out of 21 tested cytomegalovirus (CMV) specimens.

The CheckNOW™ HIV SELF TEST was also evaluated with 23 interfering substances which include medicine and blood analyte. These substances were spiked with HIV-1 Antibody positive plasma and the test results indicated these interference substances did not affect the performance of the CheckNOW™ HIV SELF TEST.

LITERATURE REFERENCES

- Blattner, W., Gallo, R., & Temin, H. HIV causes AIDS. Science. 1988; 241(4865): 515-515.
- CDC: 2008 Case Definition; Human Immunodeficiency Virus Infection.
- Delaney KP, Branson BM, Uniyal A, et al. Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. Clinical Infectious Diseases. 2011; 52(2): 257-263.
- O'Connell RJ, Merritt TM, Malia JA, et al. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. Journal of Clinical Microbiology. 2003; 41(5): 2153-2155.
- O'Connell RJ, Agan BK, Anderson SA, et al. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. Journal of Clinical Microbiology.

TECHNICAL SUPPORT

Europe & Middle East Tel. +44 161 483 9032
Africa Tel. +27 10 500 9700
Asia Pacific Tel. +61 7 3363 7100
Taiwan Tel. 00801128036
Russia & CIS Tel. +7 499 403 9512
Latin America Tel. +57 1 482 4033
Brazil Tel. 0800 111 3363

Email: EME.techsupport@abbott.com
Email: ARCIS.techsupport@abbott.com
Email: AP.techsupport@abbott.com
Email: TW.techsupport@abbott.com
Email: ARCIS.techsupport@abbott.com
Email: LA.techsupport@abbott.com
Email: LA.techsupport@abbott.com

Abbott Rapid Diagnostics Jena GmbH
Orlweg 1, D-07743 Jena, Germany
www.abbott.com/poct

©2022 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.