

Supplementary information for the manuscript titled
Association of *Wnt4* (rs121908652) and *Bmp15* (rs104894766) polymorphisms with endometriosis
risk in Iranian women

In this study, the ready-to-use GeneRuler™ 50 bp DNA Ladder molecular marker (Zymogram Gels, Thermo Fisher Scientific, USA) was used. This marker has a concentration of 0.1 µg/µL and is stored in a special marker storage buffer. This buffer contains 10 mM EDTA and 10 mM Tris-HCl, stored at -20°C. After horizontal electrophoresis on agarose gel, the marker fragments were separated from each other and positioned at 50 bp intervals (ranging from 50–1000 bp) apart (Figure S1).




Figure S1 Fragments obtained from the GeneRuler™ 50 bp DNA Ladder