

Supplementary Table 1: Details of the quantitative real-time PCR primers

Gene Name	Forward Primer	Reverse Primer
<i>LAMC2</i>	5'- ACCGAGTTCGGGATAC-3' 5'-CCAGAAAGGTTTCAGATGCC-3'	5'-TGCTCCATGTTACTGCG-3' 5'-TGGCTTCCAAGTTCAGACTC-3'
<i>Vimentin</i>	5'-CTGCCAACCGGAACAATGAC-3' 5'-TGCAGGAGGAGATGCTTCAG-3'	5'-CATTTACGCATCTGGCGTT-3' 5'-ATTCCACTTTGCGTTCAAGG-3'
<i>ZEB1</i>	5'-AGGATGACCTGCCAACAGAC-3' 5'-GCCAATAAGCAAACGATTCTG-3'	5'-CTTCAGGCCCCAGGATTTCTT-3' 5'-TTTGGCTGGATCACTTTCAAG-3'
<i>SNAIL1</i>	5'-TCTTTCCTCGTCAGGAAGCC-3'	5'-GATCTCCGGAGGTGGGATGG-3'
<i>E-Cadherin</i>	5'-TTACTGCCCCCGAGAGGAGGA-3' 5'-AGAACGCATTGCCACATACACTCTC-3'	5'-TGCAACGTCGTTACGAGTCA-3' 5'-CGGTTACCGTGATCAAAATCTCCA-3'
<i>N-Cadherin</i>	5'-CCGGTTTCATTTGAGGGCAC-3'	5'-TCCCTCAGGAACTGTCCAT-3'
<i>GAPDH</i>	5'-GGAAGGTGAAGGTCGGAGT-3'	5'-CCATGGGTGGAATCATATTGGA-3'
<i>U6</i>	5'-GTGCTCGCTTCGGCAGCACATATAC-3'	5'-AAAAATATGGAACGCTTCACGAATTTG-3'