

Supplemental Table 2 Fluorescence antibody used for flow cytometry and sorting				
Fluorochrome	Antigen	Clone	Supplier	Category
FVS780	Cell surface and intracellular amines	/	BD	discrimination of live/dead cells
FITC	CD4	RPA-T4	BioLegend	surface
PE/Cyanine7	CD45RA	HI100	BioLegend	surface
APC	CD45RO	UCHL1	BioLegend	surface
PE	SIGIRR	Clone #010	Sino Biological	surface
BUV395	TNF- α	MAB11	BD	intracellular
BV421	IL-17A	N49-653	BD	intracellular
BV510	CD45	30-F11	BD	surface
APC	CD11b	M1/70	BD	surface
FITC	Gr-1	RB6-8C5	BioLegend	surface
PE	F4/80	T45-2342	BD	surface
Alexa Fluor 647	TNF- α	MP6-XT22	BioLegend	intracellular
PE/Cyanine7	IL-17A	TC11-18H10.1	BioLegend	intracellular