

Oxidation driven reversal of PIP₂-dependent gating in GIRK2 channels

Sun-Joo Lee* and Colin G. Nichols

¹Department of Cell Biology and Physiology and the Center for Investigation of Membrane Excitability Diseases, Washington University School of Medicine, St. Louis, Missouri, USA,


To whom correspondence should be addressed.

Key Words: GIRK2, Kir3, PIP₂, oxidation, oxidative stress, reversed gating.

Correspondence:


Ph: +1 314 362 6629, FAX: +1 314 362 7463, email: sunjoo.lee@wustl.edu

Supplementary information




Supple Fig. S1. Purified GIRK2 proteins in reducing and oxidizing condition

The elution profiles of GIRK2 wild type protein (**A**), GIRK2 mutants (C65V and C221T/C321V) (**B**), and C190S mutant (**C**) off the size exclusion columns equilibrated in reducing (with 2 mM TCEP; black) or in oxidizing (without TCEP; red) condition. Coomassie stained 1D SDS-PAGE image of the each sample. The type of the used size exclusion column is noted above the profiles.





Supple Fig. S2 K⁺ specific basal activity elevated in oxidized GIRK2 proteins

Differential ACMA signal decay (ΔF) of GIRK2 proteins in the oxidizing condition, purified and incubated for 3 hours without TCEP to test basal channel activity with increasing amount of proteins and with the internal salt KCl (Left) or NaCl (Middle). Comparison of K⁺ vs Na⁺ specific basal activities (Right).



Supple Fig. S3. Purified cKir2.2 protein in reducing and oxidizing condition

A. The elution profiles of cKir2.2 off the size exclusion column equilibrated in reducing (+2 mM TCEP; black) or oxidizing (0 mM TCEP; red) condition, and the Coomassie stained 1D SDS-PAGE profiles of the protein samples.



Supple Fig. S4. cKir2.2 channel activities in the reducing and oxidizing condition
A. Differential ACMA signal decay (ΔF) by cKir2.2 purified with 2 mM TCEP in reducing condition, and activities shown under different combinations of agonists with external NMDG (blue) or NaCl (green) and with 0 (open) or 0.1 (filled) % PIP2. **B. Left:** Differential ACMA signal decay (ΔF) by cKir2.2 purified without TCEP in oxidizing condition, incubated for 3 hours either with (black) or without (red) TCEP. The buffer composition for the ACMA assay with internal NaCl and external NMDG is shown in the schematic diagram. **Right:** Activity in oxidizing condition is in red and in reducing condition in black. Paired student t-test was applied with $n = 3$.