

MITS Ambulance sterilisation

- Roof, floor and 3 walls of MITS ambulance are made air tight.
- Roof, floor, all walls and door are made smooth to prevent accumulation of dust.
- The MITS ambulance is cleaned and washed thoroughly from inside and outside.
- There is regular cleaning of all inner surfaces of ambulance with Bacillol 25 with spray, at least once a week and after MITS.
- MITS ambulance is always ready by fogging with Glutaral (15.2gm) and Dimethanol (19.7gm) (after MITS and before MITS) to maintain asepsis .
- Fogging with H₂O₂ , Silver nitrate for one hour is done during emergency or Air sterilisers are used for instant sterilisation before MITS procedure.
- Air conditioner is installed in MITS ambulance.