

Table 1. The primers sequences for ChIP assay.

Name	Sequences (5' - 3')	Primer Length (bp)
ChIP-B7BS-F1	ACACCTGGCTCATTTTTCTT	20
ChIP-B7BS-R1	AAAATTAGCTGGGTATGGTG	20
ChIP-B7BS-F2	TGTATTTT TAGTAGAGTGGG	20
ChIP-B7BS-R2	TTTGGAACCTGAATGAATT	20
ChIP-B7BS-F3	CAACCAACTGGCAATAAACTC	21
ChIP-B7BS-R3	GCATGTGTTAAATAATGTTGT	21
ChIP-B7BS-F4	AATAAAATGGGAGAGGGGGTG	21
ChIP-B7BS-R4	AGGACTACAGACTGCCTCTAG	21