

Table S2 Primers used in quantitative real-time PCR (qRT-PCR) assay

Primer	Sequence (5'-3')
BMI-1-Forward	CCACAGTTTCCTCACATTTCC
BMI-1-Reverse	AAGTTGCTGATGACCCATTAC
CD44-Forward	CTGCCGCTTTGCAGGTGTA
CD44- Reverse	CATTGTGGGCAAGGTGCTATT
CD133-Forward	TTCTTGACCGACTGAGACCCA
CD133- Reverse	TCATGTTCTCCAACGCCTCTT
SOX2-Forward	TGGACAGTTACGCGCACAT
SOX2-Reverse	CGAGTAGGACATGCTGTAGGT
OCT4-Forward	CTTGAATCCCGAATGGAAAGGG
OCT4-Reverse	GTGTATATCCCAGGGTGATCCTC
MYCN-Forward	ACCCGGACGAAGATGACTTCT
MYCN-Reverse	CAGCTCGTTCTCAAGCAGCAT
GAPDH-Forward	CAACGGATTTGGTCGTATTGG
GAPDH- Reverse	TGATGGCAACAATATCCACTTTACC