

Supp. Table 2. Antibodies used in the study

Name	Dilution	Manufacturer
Cone arrestin	1:5000	Millipore, Billerica, MA
GFAP	1:1000	Sigma, St. Louis, MO
Ki67	1:1000	ThermoScientific
MAB1281	1:300	Millipore, Billerica, MA
Nestin	1:10000	Millipore, Billerica, MA
PKCa	1:8000	Sigma, St. Louis, MO
Recoverin	1:2000	Millipore, Billerica, MA
Stem 121	1:300	Stem Cells, Newark, CA
Stem 123	1:3000	Stem Cells, Newark, CA
Synaptophysin	1:1000	Millipore, Billerica, MA