

Article title: Determination of cortisol cut-off limits and steroid dynamics in ACTH stimulation (Synacthen®) test: A comparative analysis using Roche Elecsys Cortisol II immunoassay and LC-MS/MS

Journal name: Endocrine (Springer)

Author names: Sema Okutan^{1,2}, Nanna Thurmann Jørgensen^{1,2}, Lars Engers Pedersen³, Stina Willemoes Borresen¹, Linda Hilsted⁴, Lennart Friis Hansen^{3,5}, Ulla Feldt-Rasmussen^{1,2}, Marianne Klose¹

Affiliations:

¹Department of Endocrinology and Metabolism, Copenhagen University Hospital, Rigshospitalet, Copenhagen, Denmark

²Department of Clinical Medicine, Faculty of Health and Medical Sciences, Copenhagen University, Copenhagen, Denmark

³Department of Clinical Biochemistry, Næstved, Slagelse and Ringsted Hospitals, Slagelse, Denmark

⁴Department of Clinical Biochemistry, Copenhagen University Hospital, Rigshospitalet, Copenhagen, Denmark

⁵Department of Clinical Biochemistry, Copenhagen University Hospital, Bispebjerg Hospital, Copenhagen, Denmark

Corresponding author's e-mail address: marianne.christina.klose.01@regionh.dk

Online Resource 1. The LC-MS/MS method.

	Ion	Daughter	Dwell(s)	Cone(V)	Coll(eV)
Cortisol 1	363.3	97.00	0.014	44.00	30.00
Cortisol 2	363.3	121.0	0.014	44.00	32.00
Internal standard	367.3	121.0	0.014	44.00	24.00