

Supplemental Table 2. The primary antibodies and primer sequences information of multiple genes for Western blot and PCR experiments.

Genes	Primary antibodies	Primer sequences	
		Forward primer	Reverse primer
RPL17	dilutions with 1: 1,500, Abcam, USA.	CGTTAATCGCTATCGGAATC	TTGCTGAACATCAGTATTA
RPS9	dilutions with 1: 1,000, Abcam, USA.	TTGGCCTCGCCTTTGCCCGA	TGCACACGCCCTGGAATGC
P53	dilutions with 1: 4,000, Proteintech, China.	CTCCGGCATCTTATCCGAGTGG	TATTAGGTGGTACAGTCAGAGC
Bax	dilutions with 1: 1,000, Proteintech, China.	ATCGGATGCGTCCACCAAATTG	ATTGCCCACGGCGGCACGTTTC
BCL2	dilutions with 1: 2,000, Proteintech, China.	GTTACCCTGTGGATGACTCCAT	CCCTAGAGAAATCAAACAGCGTA
Cyclin D	dilutions with 1: 1,000, Proteintech, China.	CGTAAACCGACAACCTCCAAAAT	AATGGCATT TTTGGAGAGGACGTT