

Resident Name:

Date:

Initial setup: Cardiac preset, minimum depth, gain very high									Score	
Time limit: 8 minutes									Total	Max
<b>Probe selection</b> - "Please choose the most appropriate probe for a focused lung and abdominal exam")										1 point
<ul style="list-style-type: none"> <li>1 point for phased array probe</li> </ul>										
<b>Preset selection</b> - "Please select the exam preset on the machine for a focused lung exam"										1 point
<ul style="list-style-type: none"> <li>1 point for lung (or abdomen if preference / no lung preset)</li> </ul>										
	Structure		Gain		Depth		Probe marker			
	Score	Max	Score	Max	Score	Max	Score	Max		
<b>Lung sliding</b>		2 points • Ribs • Pleura		2 points		1 point		1 point		6 points
<b>Lung A-lines</b>		1 point • Min 2 a-lines				1 point				2 points

<b>Preset selection</b> - "Please select the exam preset on the machine for a focused abdominal exam"										1 point
<ul style="list-style-type: none"> <li>1 point for abdomen preset (or if no change required)</li> </ul>										
	Structure		Gain		Depth		Probe marker			
	Score	Max	Score	Max	Score	Max	Score	Max		
<b>RUQ</b>		3 points • Diaphragm • Liver • kidney		2 points		1 point		1 point		7 points
<b>LUQ</b>		3 points • Diaphragm • Spleen • kidney				1 point		1 point		5 points
<b>Kidney longitudinal</b>		1 point • Full view pole to pole				1 point		1 point		3 points
<b>Transverse bladder</b>		1 point • Max dimension		2 points		1 point		1 point		5 points
<b>Sagittal bladder</b>		1 point • Max dimension				1 point		1 point		3 points

<b>Total Score: / 34</b>
--------------------------

Gain score:

- 1 point - Some adjustment in appropriate direction, but inadequate/suboptimal
- 2 points - Gain optimized (or if no adjustment necessary)

Depth score:

- 1 point - Structures of interest take up >2/3 of screen (or if no adjustment necessary)

Probe marker score:

- 1 point - Probe marker should be towards head unless specifically indicated by view