




Supplementary Figure S1. Identification of alveolar and splenic macrophages. Gating strategy to identify macrophages in the lung (**a**), and spleen (**b**). Representative plots are shown. The results are representative of three independent experiments. Mφ, macrophage.

Supplementary Table S1 Antibodies used for flow cytometry

Antigen detected	Clone	Manufacturer
Annexin V		BioLegend
B220 (CD45R)	RA3-6B2	BioLegend
CD11b	M1/70	BioLegend
CD11c	N418	BioLegend
CD14	Sa14-2	BioLegend
CD45	30-F11	BioLegend
CD45.1	A20	BioLegend
CD45.2	104	BioLegend
CD64	X54-5/7.1	BioLegend
CD206	C068C2	BioLegend
CX3CR1	2A9-1	BioLegend
F4/80	BM8	BioLegend
MHCII (I-A/I-E)	M5/114.15.2	BioLegend
MerTK	Mer	BioLegend
Siglec-F	E50-2440	BD Pharmingen

Supplementary Table S2 Primers used for quantitative RT-PCR

Primer	Sequence
<i>b-actin</i> Forward	5'-CGTTGACATCCGTAAAGACCTC-3'
<i>b-actin</i> Reverse	5'-AGCCACCGATCCACACAGA-3'
<i>Csf1r</i> Forward	5'-AGGGCCATATACAGGTACACAT-3'
<i>Csf1r</i> Reverse	5'-CACAGGCATCCATGTAACAC-3'
<i>Csf2rb</i> Forward	5'-GGGCAGGAACACAGGACTTCAGGAC-3'
<i>Csf2rb</i> Reverse	5'-CCAGGCCTCTAGCTACCTTGACAGG-3'
<i>Csf1</i> Forward	5'-TACAAGTGGAAGTGGAGGAGCCAT-3'
<i>Csf1</i> Reverse	5'-AGTCCTGTGTGCCCAGCATAGAAT-3'
<i>Csf2</i> Forward	5'-GGCTAAGGTCCTGAGGAGGAT-3'
<i>Csf2</i> Reverse	5'-ACCTCTTCATTCAACGTGACAGG-3'