

# DICOM Validation Report – 2025-09-13 19:13

## File Summary

File	Errors	Warnings	Private	Total
All files	23	10	0	33
DCM_test_CT_각시_slice1177.dcm	23	10	0	33

## Detailed Messages

File	Severity	Tag	Message
DCM_test_CT_각시_slice1177.dcm	ERR		Error - MediaStorageSOPInstanceUID different from SOPInstanceUID
DCM_test_CT_각시_slice1177.dcm	ERR		Error - PixelAspectRatio may not be present when it has a ratio of 1:1 - values are 1\1
DCM_test_CT_각시_slice1177.dcm	ERR	0018,1150	Error - Value invalid for this VR - (0x0018,0x1150) IS Exposure Time IS [1] = <0.000000> - Character invalid for this VR = '.' (0x2e)
DCM_test_CT_각시_slice1177.dcm	ERR	0018,1151	Error - Value invalid for this VR - (0x0018,0x1151) IS X-Ray Tube Current IS [1] = <0.000000> - Character invalid for this VR = '.' (0x2e)
DCM_test_CT_각시_slice1177.dcm	ERR	0018,1152	Error - Value invalid for this VR - (0x0018,0x1152) IS Exposure IS [1] = <0.000000> - Character invalid for this VR = '.' (0x2e)
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Dicom dataset contains invalid data values for Value Representations
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 2 Required Element= <PatientID> Module= <Patient>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 2 Required Element= <PatientBirthDate> Module= <Patient>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 2 Required Element= <PatientSex> Module= <Patient>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 2 Required Element= <ReferringPhysicianName> Module= <GeneralStudy>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 2 Required Element= <StudyID> Module= <GeneralStudy>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 2 Required Element= <AccessionNumber> Module= <GeneralStudy>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 2 Required Element= <SeriesNumber> Module= <GeneralSeries>

File	Severity	Tag	Message
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 2C Conditional Element= <Laterality> Module= <GeneralSeries>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 2C Conditional Element= <PatientPosition> Module= <GeneralSeries>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 1 Required Element= <FrameOfReferenceUID> Module= <FrameOfReference>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 2 Required Element= <PositionReferenceIndicator> Module= <FrameOfReference>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Shall not be present when PixelSpacing is present - attribute <PixelAspectRatio>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 1 Required Element= <ImageType> Module= <CTImage>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - A value is required for value 3 in CT Images - attribute <ImageType>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 1 Required Element= <RescaleIntercept> Module= <CTImage>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 1 Required Element= <RescaleSlope> Module= <CTImage>
DCM_test_CT_각시_slice1177.dcm	ERR		Error - Missing attribute Type 2 Required Element= <AcquisitionNumber> Module= <CTImage>
DCM_test_CT_각시_slice1177.dcm	WARN		Warning - Missing attribute or value that would be needed to build DICOMDIR - Patient ID
DCM_test_CT_각시_slice1177.dcm	WARN		Warning - Missing attribute or value that would be needed to build DICOMDIR - Study ID
DCM_test_CT_각시_slice1177.dcm	WARN		Warning - Missing attribute or value that would be needed to build DICOMDIR - Series Number
DCM_test_CT_각시_slice1177.dcm	WARN	0010,0010	Warning - Value dubious for this VR - (0x0010,0x0010) PN Patient's Name PN [1] = <SfxImage2D> - Retired Person Name form
DCM_test_CT_각시_slice1177.dcm	WARN		Warning - Value is zero for value 1 of attribute <KVP>
DCM_test_CT_각시_slice1177.dcm	WARN		Warning - Value is zero for value 1 of attribute <Distance Source to Detector>
DCM_test_CT_각시_slice1177.dcm	WARN		Warning - Value is zero for value 1 of attribute <Distance Source to Patient>
DCM_test_CT_각시_slice1177.dcm	WARN		Warning - Value is zero for value 1 of attribute <Exposure Time>

File	Severity	Tag	Message
DCM_test_CT_각시_slice1177.dcm	WARN		Warning - Value is zero for value 1 of attribute <X-Ray Tube Current>
DCM_test_CT_각시_slice1177.dcm	WARN		Warning - Value is zero for value 1 of attribute <Exposure>

## Full dciodvfy Dumps

▶ DCM\_test\_CT\_각시\_slice1177.dcm

## DICOM Tag Logs (from dataset)

▶ DCM\_test\_CT\_각시\_slice1177.dcm

Raw JSONL: C:\Users\USER\Downloads\20250912\_CT 규격dcm과 기존 dcm 제작\_태그 샘플 작성  
\output\_20250913\_1913\validation\_20250913\_191355\dciodvfy.jsonl