

Supplementary Table 1 - List of antibodies used in this study

Panel	Target	Clone/Reference	Supplier	Isotype	Dye
Myeloid and T-cell populations	MHCII	ILA21	In house	IgG2a	BV605
	CD3	MM1A	Clinisciences	IgG1	BV711
	G1	CH138A	Kingfisher	IgGM	PE-Cy7
	CD14	Tuk4	Biorad	IgG2a	RPE
	Viability	65-0862-14	Thermo Fisher Scientific		eFluor 450
T-cell populations	CD3	MM1A	Clinisciences	IgG1	BV711
	CD4	CC8	Biorad	IgG2a	A647
	CD8	CC63	Biorad	IgG2a	BV605
	TCR- $\gamma\delta$	GB21A	Clinisciences	IgG2b	BV510
	IL-17	eBio64DEC17	Fisher scientific	IgG1	PE-Cy7
	IFN γ	CC302	Biorad	IgG1	A488
	Viability	65-0862-14	Thermo Fisher Scientific		eFluor 450
Cell sorting for scRNA analysis	G1	CH138A	Kingfisher	IgGM	
	CD45	MCA2220PE	Biorad	IgG1	RPE
	Viability	65-0862-14	Thermo Fisher Scientific		eFluor 450