

Supplementary Information

Dry-printed electrodes for transcutaneous electrical stimulation of innervated muscles: Towards wearable and closed-loop stimulation

Rawan Ibrahim^{1,2}, Ieva Vėbraitė^{1,2}, Daniele Re^{1,2}, Paul F. Funk^{3,4}, Gerd Fabian Volk^{3,4}, Orlando Guntinas-Lichius^{3,4}, and Yael Hanein^{1,2,5,6,*}

¹School of Electrical and Computer Engineering, Tel Aviv University, Tel Aviv, Israel

²Tel Aviv University Center for Nanoscience and Nanotechnology, Tel Aviv, Israel

³Facial-Nerve-Center Jena, Jena University Hospital, Jena, Germany

⁴Center for Rare Diseases, Jena University Hospital, Jena, Germany

⁵Sagol School of Neuroscience, Tel Aviv University, Tel Aviv, Israel

⁶X-trodes, Herzlia, Israel

*yaelha@tauex.tau.ac.il

Supplementary figures




Figure S1. Setup for transcutaneous electrical stimulation using dry-printed electrodes. The system includes (1) a stimulus generator (STG5, Multi Channel Systems), (2) BNC-to-1 mm adaptors, (3) a custom PCB connector, and (4) a soft printed electrode array. This configuration ensures stable electrical contact while maintaining flexibility and skin conformity.

Supplementary Videos

Video 1: Demonstration of accessory nerve stimulation using dry and gel electrodes.

Video 2: Demonstration of closed-loop facial stimulation.