

Examples of original Western blot for three repeats

1.The target protein band are circled with blue dotted boxes.

2.To conduct the experiments better, we cut the entire membrane into several parts according to the molecular protein and then spliced them together during imaging.

3.Each slide shows the results of three independent repeats of the original protein bands of a certain target protein weight marked)




Figure 4 A




Figure 4 A




Figure 5 A




Figure 5 A




Figure 5 E




Figure 5 E




Figure 7 A




Figure 7A




Figure 8 A




Figure 8 A




Figure 8 D




Figure 8 D




Figure 8 D